

A-Line™ 自動調芯レンズマウント



複数のレンズを素早く簡単にマウント&アライメント

各レンズマウントの 6mm アライメント穴にガイドロッドを通して連結するだけで直列アライメントが完了！

- **多彩なマウント：** オプティクス径 0.5～3 インチ、C-Mount/RMS スレッド/ミラー用二重穴/30mm レンズ用
- **効率的なデザイン：** スリムエッジ、大開口、スレッドポスト・マウント
- **簡単セットアップ：** 自己整列、素早いフォーカシング、グループでマウント可能、近接設置
- **高精度ガイドロッド：** 6mm 径 4-40 スレッド、1～24 インチ長、目盛り付きまたは回転方向を揃えるキー加工



株式会社日本レーザー

東京本社 〒169-0051 東京都新宿区西早稲田 2-14-1 TEL 03-5285-0853
大阪支店 〒533-0033 大阪市東淀川区東中島 1-20-12 TEL 06-6323-7286
名古屋支店 〒460-0003 名古屋市中区錦 3-1-30 TEL 052-205-9711
e-mail: newport@japanlaser.jp / <http://www.newport-japan.jp>

A-Line 固定オプティクスマウント特性表

Model	Optic Diameter [in. (mm)]	Clear Aperture [in. (mm)]	Lens Thread:	Optical Axis Height [in. (mm)]	Max. Lens Thk.(mm)	Post Thd.
M-LH-0.5A	0.5 (12.7)	0.43 (10.9)	0.535-40	1.0 (25.4)	6	M4
M-LH-1A	1.0 (25.4)	0.90 (22.9)	1.035-40	1.25 (31.8)	7	M4
M-LH-1T	1.0 (25.4)	0.90 (22.9)	1.035-40	1.25 (31.8)	3	M4
M-LH-30M	1.2 (30.0)	1.10 (27.9)	M30.5 X 0.5	1.25 (31.8)	7	M4
M-LH-1.5A	1.5 (38.1)	1.40 (35.6)	1.535-40	1.5 (38.1)	9	M4
M-LH-2A	2.0 (50.8)	1.87 (47.5)	2.035-40	1.75 (44.5)	9	M4
M-LH-3A	3.0 (76.2)	2.80 (71.1)	3.035-40	2.25 (57.2)	9	M4

A-Line マウントアダプター一覧

Model	Description
LH-1PM	A- LINE Fixed Cube Mount, 1.0 in. (25.4 mm) , M4 counter bore
M-LH-1DB	A-LINE Fixed Double-Bore Mount, 1.0 in. (25.4 mm) Diameter, M4 Thd.
M-LH-CMT	A-LINE Fixed Lens Mount, C-Mount Camera Lens Thread (1.00-32), M4 Thd.
M-LH-K1	A-LINE Insertable Link Block, 1/2-in to 1-in, M4 Thd.
M-LH-K2	A-LINE Insertable Link Block, 1-in to 2-in, M4 Thd.
M-LH-RMS	A-LINE Fixed Microscope Objective (0.800-36 Thd.) Holder M4 Thd.



GR シリーズ ガイドロocker一覧

Model	Description
LH-GR-0.5	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 0.5 in. long
LH-GR-1	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 1 in. long
LH-GR-10	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 10 in. long
LH-GR-12	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 12 in. long
LH-GR-2	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 2 in. long
LH-GR-3	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 3 in. long
LH-GR-4	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 4 in. long
LH-GR-6	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 6 in. long
LH-GR-8	A-LINE D-Keyed Guide Rod, 6 mm, 4-40 Thd, 8 in. long
OC1-GR-1	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 1 in. long
OC1-GR-10	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 10 in. long
OC1-GR-12	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 12 in. long
OC1-GR-18	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 18 in. long
OC1-GR-2	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 2 in. long
OC1-GR-24	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 24 in. long
OC1-GR-3	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 3 in. long
OC1-GR-4	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 4 in. long
OC1-GR-6	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 6 in. long
OC1-GR-8	1 in. OpticsCage+ Guide Rod, 6 mm, 4-40 Thd, 8 in. long
OC1-GR-4S	1 in. OpticsCage+ Guide Rods, 6 mm, 4-40 Thd, 4 in. long, 1 mm laser marked
OC1-GR-6S	1 in. OpticsCage+ Guide Rods, 6 mm, 4-40 Thd, 6 in. long, 1 mm laser marked



株式会社日本レーザー

東京本社 〒169-0051 東京都新宿区西早稲田 2-14-1 TEL 03-5285-0853
 大阪支店 〒533-0033 大阪市東淀川区東中島 1-20-12 TEL 06-6323-7286
 名古屋支店 〒460-0003 名古屋市中区錦 3-1-30 TEL 052-205-9711
 e-mail: newport@japanlaser.jp / <http://www.newport-japan.jp>